





TBC1D31 Antibody

Product Code	CSB-PA822197LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96DN5
Immunogen	Recombinant Human TBC1 domain family member 31 protein (662-1001AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TBC1 domain family member 31 (WD repeat-containing protein 67), TBC1D31, WDR67
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TBC1D31

Image



Immunohistochemistry analysis of human testis tissue using CSB-PA822197LA01HU at dilution of 1:100



Immunocytochemistry analysis of human glioma using CSB-PA822197LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.