




NECTIN2 Antibody

Product Code	CSB-PA835690LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92692
Immunogen	Recombinant Human Nectin-2 protein (32-360AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nectin-2 (Herpes virus entry mediator B) (Herpesvirus entry mediator B) (HveB) (Nectin cell adhesion molecule 2) (Poliovirus receptor-related protein 2) (CD antigen CD112), NECTIN2, HVEB PRR2 PVRL2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NECTIN2

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, A549 whole cell lysate, HeLa whole cell lysate

All lanes: NECTIN2 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 58, 52 kDa

Observed band size: 58 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.