



TSR2 Antibody, HRP conjugated

Product Code	CSB-PA836166HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969E8
Immunogen	Recombinant Human Pre-rRNA-processing protein TSR2 homolog protein (1-191AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pre-rRNA-processing protein TSR2 homolog, TSR2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TSR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.