



AKT1S1 Antibody, Biotin conjugated

Product Code	CSB-PA836185LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96B36
Immunogen	Recombinant Human Proline-rich AKT1 substrate 1 protein (52-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Proline-rich AKT1 substrate 1 (40 kDa proline-rich AKT substrate), AKT1S1, PRAS40
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	AKT1S1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.