



VPS39 Antibody, Biotin conjugated

Product Code	CSB-PA836235LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JC1
Immunogen	Recombinant Human Vam6/Vps39-like protein (263-490AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vam6/Vps39-like protein (TRAP1-like protein) (hVam6p), VPS39, KIAA0770 TLP VAM6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS39