



VPS39 Antibody, Biotin conjugated

Product Code	CSB-PA836235LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JC1
Immunogen	Recombinant Human Vam6/Vps39-like protein (263-490AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vam6/Vps39-like protein (TRAP1-like protein) (hVam6p), VPS39, KIAA0770 TLP VAM6
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VPS39
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.