





ELMO2 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA836237LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96JJ3 |
| Immunogen | Recombinant Human Engulfment and cell motility protein 2 protein (304-371AA) |
| Raised In | Rabbit |
| Species Reactivity | Human, Mouse |
| Tested Applications | ELISA, WB, IHC, IF; Recommended dilution: WB:1:200-1:5000, IHC:1:20-1:200, IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Engulfment and cell motility protein 2 (Protein ced-12 homolog A) (hCed-12A), ELMO2, CED12A KIAA1834 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ELMO2 |


Image



Immunohistochemistry of paraffin-embedded human placenta tissue using CSB-PA836237LA01HU at dilution of 1:100



Immunofluorescent analysis of A549 cells using CSB-PA836237LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)



Western Blot

Positive WB detected in: Mouse heart tissue

All lanes: ELMO2 antibody at 3.4μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 83, 73 kDa

Observed band size: 83 kDa