



NOA1 Antibody, Biotin conjugated

Product Code	CSB-PA836715LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NC60
Immunogen	Recombinant Human Nitric oxide-associated protein 1 protein (42-174AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nitric oxide-associated protein 1, NOA1, C4orf14
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NOA1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.