



Cd207 Antibody

Product Code	CSB-PA837695LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VBX4
Immunogen	Recombinant Mouse C-type lectin domain family 4 member K protein (63-331AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-type lectin domain family 4 member K (Langerin) (CD antigen CD207), Cd207, Clec4k
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Cd207