



GUCD1 Antibody, Biotin conjugated

Product Code	CSB-PA839384LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96NT3
Immunogen	Recombinant Human Protein GUCD1 protein (1-136AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein GUCD1 (Guanylyl cyclase domain-containing protein 1) (Protein LLN4), GUCD1, C22orf13 LLN4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GUCD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.