



DME Antibody, Biotin conjugated

Product Code	CSB-PA839701LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8LK56
Immunogen	Recombinant Arabidopsis thaliana Transcriptional activator DEMETER protein (1-320AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional activator DEMETER (EC 3.2.2.-) (DNA glycosylase-related protein DME), DME
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	DME
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.