



MYRIP Antibody, HRP conjugated

Product Code	CSB-PA839876LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NFW9
Immunogen	Recombinant Human Rab effector MyRIP protein (387-689AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab effector MyRIP (Exophilin-8) (Myosin-VIIa- and Rab-interacting protein) (Synaptotagmin-like protein lacking C2 domains C) (SlcC2-c) (Slp homolog lacking C2 domains c), MYRIP, SLAC2C
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	MYRIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.