




TAS1R2 Antibody

Product Code	CSB-PA840577LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE23
Immunogen	Recombinant Human Taste receptor type 1 member 2 protein (183-362AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IF; Recommended dilution: IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 1 member 2 (G-protein coupled receptor 71) (Sweet taste receptor T1R2), TAS1R2, GPR71 T1R2 TR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS1R2

Image



Immunofluorescent analysis of A549 cells using CSB-PA840577LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)