



TAS1R2 Antibody, HRP conjugated

Product Code	CSB-PA840577LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE23
Immunogen	Recombinant Human Taste receptor type 1 member 2 protein (183-362AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 1 member 2 (G-protein coupled receptor 71) (Sweet taste receptor T1R2), TAS1R2, GPR71 T1R2 TR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS1R2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.