



TAS1R2 Antibody, Biotin conjugated

Product Code	CSB-PA840577LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TE23
Immunogen	Recombinant Human Taste receptor type 1 member 2 protein (183-362AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Taste receptor type 1 member 2 (G-protein coupled receptor 71) (Sweet taste receptor T1R2), TAS1R2, GPR71 T1R2 TR2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	TAS1R2