



NP Antibody, Biotin conjugated

Product Code	CSB-PA841983HD01RCI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q91DE1
Immunogen	Recombinant Reston ebolavirus Nucleoprotein protein (414-739AA)
Raised In	Rabbit
Species Reactivity	Reston ebolavirus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nucleoprotein (Nucleocapsid protein) (Protein N), NP
Immunogen Species	Reston ebolavirus
Research Area	Others
Target Names	NP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.