




PIGK Antibody

Product Code	CSB-PA842143LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92643
Immunogen	Recombinant Human GPI-anchor transamidase protein (91-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI-anchor transamidase (GPI transamidase) (EC 3.-.-.-) (GPI8 homolog) (hGPI8) (Phosphatidylinositol-glycan biosynthesis class K protein) (PIG-K), PIGK, GPI8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIGK

Image



Western Blot

Positive WB detected in: HepG2 whole cell lysate, A549 whole cell lysate


All lanes: PIGK antibody at 2.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 46, 37 kDa

Observed band size: 46 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA842143LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.