



PIGK Antibody, Biotin conjugated

Product Code	CSB-PA842143LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92643
Immunogen	Recombinant Human GPI-anchor transamidase protein (91-293AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI-anchor transamidase (GPI transamidase) (EC 3.-.-.-) (GPI8 homolog) (hGPI8) (Phosphatidylinositol-glycan biosynthesis class K protein) (PIG-K), PIGK, GPI8
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIGK
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.