





USP6NL Antibody

Product Code	CSB-PA842149LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92738
Immunogen	Recombinant Human USP6 N-terminal-like protein (387-693AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	USP6 N-terminal-like protein (Related to the N-terminus of tre) (RN-tre), USP6NL, KIAA0019
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	USP6NL

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA842149LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA842149LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.