



LHFPL5 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA844703LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8TAF8 |
| Immunogen | Recombinant Human LHFPL tetraspan subfamily member 5 protein (46-99AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | LHFPL tetraspan subfamily member 5 protein (Lipoma HMGIC fusion partner-like 5 protein) (Tetraspan membrane protein of hair cell stereocilia), LHFPL5, TMHS |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | LHFPL5 |