



LHFPL5 Antibody, FITC conjugated

Product Code	CSB-PA844703LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAF8
Immunogen	Recombinant Human LHFPL tetraspan subfamily member 5 protein (46-99AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LHFPL tetraspan subfamily member 5 protein (Lipoma HMGIC fusion partner-like 5 protein) (Tetraspan membrane protein of hair cell stereocilia), LHFPL5, TMHS
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	LHFPL5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.