



GABPB2 Antibody, HRP conjugated

Product Code	CSB-PA844704LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAK5
Immunogen	Recombinant Human GA-binding protein subunit beta-2 protein (163-440AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GA-binding protein subunit beta-2 (GABP subunit beta-2) (GABPB-2), GABPB2
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GABPB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.