




# Loxl4 Antibody

<b>Product Code</b>	CSB-PA846069LA01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q924C6
<b>Immunogen</b>	Recombinant Mouse Lysyl oxidase homolog 4 protein (26-757AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse, Rat
<b>Tested Applications</b>	ELISA, WB; Recommended dilution: WB:1:500-1:5000
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Lysyl oxidase homolog 4 (EC 1.4.3.-) (Lysyl oxidase-like protein 4) (Lysyl oxidase-related protein C), Loxl4, Loxc
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	Loxl4

## Image



Western Blot  
 Positive WB detected in: Rat brain tissue  
 All lanes: Loxl4 antibody at 4µg/ml  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/50000 dilution  
 Predicted band size: 85 kDa  
 Observed band size: 85 kDa

## Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.