



Loxl4 Antibody, HRP conjugated

Product Code	CSB-PA846069LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q924C6
Immunogen	Recombinant Mouse Lysyl oxidase homolog 4 protein (26-757AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lysyl oxidase homolog 4 (EC 1.4.3.-) (Lysyl oxidase-like protein 4) (Lysyl oxidase-related protein C), Loxl4, Loxc
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Loxl4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.