





PREX1 Antibody

Product Code	CSB-PA848403LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCU6
Immunogen	Recombinant Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 1 protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 1 protein (P-Rex1) (PtdIns(3,4,5)-dependent Rac exchanger 1), PREX1, KIAA1415
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PREX1


Image



Immunohistochemistry analysis of human colon cancer using CSB-PA848403LA01HU at dilution of 1:100



Immunocytochemistry analysis of human lung tissue using CSB-PA848403LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA848403LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.