



PREX1 Antibody, HRP conjugated

Product Code	CSB-PA848403LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCU6
Immunogen	Recombinant Human Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 1 protein (1-44AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 3,4,5-trisphosphate-dependent Rac exchanger 1 protein (P-Rex1) (PtdIns(3,4,5)-dependent Rac exchanger 1), PREX1, KIAA1415
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PREX1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.