



ascl1a Antibody

Product Code	CSB-PA852654LA01DIL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q90259
Immunogen	Recombinant Zebrafish Achaete-scute homolog 1a protein (1-196AA)
Raised In	Rabbit
Species Reactivity	Zebrafish
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Achaete-scute homolog 1a (Zash-1a) (Pituitary-absent protein), ascl1a, ash ash1a asha pia
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ascl1a
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.