



WASF1 Antibody, HRP conjugated

Product Code	CSB-PA852883LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92558
Immunogen	Recombinant Human Wiskott-Aldrich syndrome protein family member 1 protein (381-553AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wiskott-Aldrich syndrome protein family member 1 (WASP family protein member 1) (Protein WAVE-1) (Verprolin homology domain-containing protein 1), WASF1, KIAA0269 SCAR1 WAVE1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WASF1