



WASF1 Antibody, HRP conjugated

Product Code	CSB-PA852883LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q92558
Immunogen	Recombinant Human Wiskott-Aldrich syndrome protein family member 1 protein (381-553AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Wiskott-Aldrich syndrome protein family member 1 (WASP family protein member 1) (Protein WAVE-1) (Verprolin homology domain-containing protein 1), WASF1, KIAA0269 SCAR1 WAVE1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WASF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.