



PIK3C3 Antibody, HRP conjugated

Product Code	CSB-PA854108LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEB9
Immunogen	Recombinant Human Phosphatidylinositol 3-kinase catalytic subunit type 3 protein (700-850AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositol 3-kinase catalytic subunit type 3 (PI3-kinase type 3) (PI3K type 3) (PtdIns-3-kinase type 3) (EC 2.7.1.137) (Phosphatidylinositol 3-kinase p100 subunit) (Phosphoinositide-3-kinase class 3) (hVps34), PIK3C3, VPS34
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIK3C3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.