




Metrnl Antibody

Product Code	CSB-PA855340LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8VE43
Immunogen	Recombinant Mouse Meteorin-like protein (46-311AA)
Raised In	Rabbit
Species Reactivity	Mouse, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:200-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Meteorin-like protein (Subfatin), Metrnl
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Metrnl

Image



Western Blot

Positive WB detected in: Rat brain tissue

All lanes: Metrnl antibody at 3μg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 35, 36 kDa

Observed band size: 35 kDa