



EFHD2 Antibody, Biotin conjugated

Product Code	CSB-PA856907LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C19
Immunogen	Recombinant Human EF-hand domain-containing protein D2 protein (2-240AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EF-hand domain-containing protein D2 (Swiprosin-1), EFHD2, SWS1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EFHD2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.