



PDLIM2 Antibody, Biotin conjugated

Product Code	CSB-PA856962LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JY6
Immunogen	Recombinant Human PDZ and LIM domain protein 2 protein (161-284AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ and LIM domain protein 2 (PDZ-LIM protein mystique), PDLIM2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDLIM2