





OSMR Antibody

Product Code	CSB-PA857869LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99650
Immunogen	Recombinant Human Oncostatin-M-specific receptor subunit beta protein (120-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oncostatin-M-specific receptor subunit beta (Interleukin-31 receptor subunit beta) (IL-31 receptor subunit beta) (IL-31R subunit beta) (IL-31R-beta) (IL-31RB), OSMR, OSMRB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	OSMR

Image



Immunohistochemistry of paraffin-embedded human skeletal muscle tissue using CSB-PA857869LA01HU at dilution of 1:100



Immunofluorescent analysis of MCF-7 cells using CSB-PA857869LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.