



OSMR Antibody, HRP conjugated

Product Code	CSB-PA857869LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99650
Immunogen	Recombinant Human Oncostatin-M-specific receptor subunit beta protein (120-193AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Oncostatin-M-specific receptor subunit beta (Interleukin-31 receptor subunit beta) (IL-31 receptor subunit beta) (IL-31R subunit beta) (IL-31R-beta) (IL-31RB), OSMR, OSMRB
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	OSMR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.