



GPER1 Antibody, Biotin conjugated

Product Code	CSB-PA859933LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99527
Immunogen	Recombinant Human G-protein coupled estrogen receptor 1 protein (1-62AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G-protein coupled estrogen receptor 1 (Chemoattractant receptor-like 2) (Flow-induced endothelial G-protein coupled receptor 1) (FEG-1) (G protein-coupled estrogen receptor 1) (G-protein coupled receptor 30) (GPCR-Br) (IL8-related receptor DRY12) (Lymphocyte-derived G-protein coupled receptor) (LYGPR) (Membrane estrogen receptor) (mER), GPER1, CEPR CMKRL2 DRY12 GPER GPR30
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	GPER1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.