



TRIM54 Antibody, Biotin conjugated

Product Code	CSB-PA861180LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYV2
Immunogen	Recombinant Human Tripartite motif-containing protein 54 protein (60-358AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tripartite motif-containing protein 54 (Muscle-specific RING finger protein) (MuRF) (Muscle-specific RING finger protein 3) (MuRF-3) (MuRF3) (RING finger protein 30), TRIM54, MURF MURF3 RNF30
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TRIM54
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.