



PHAX Antibody, HRP conjugated

Product Code	CSB-PA863999LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H814
Immunogen	Recombinant Human Phosphorylated adapter RNA export protein (6-243AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphorylated adapter RNA export protein (RNA U small nuclear RNA export adapter protein), PHAX, RNUXA
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PHAX
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.