



MED4 Antibody, FITC conjugated

Product Code	CSB-PA865098LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPJ6
Immunogen	Recombinant Human Mediator of RNA polymerase II transcription subunit 4 protein (167-270AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mediator of RNA polymerase II transcription subunit 4 (Activator-recruited cofactor 36 kDa component) (ARC36) (Mediator complex subunit 4) (TRAP/SMCC/PC2 subunit p36 subunit) (Vitamin D3 receptor-interacting protein complex 36 kDa component) (DRIP36), MED4, ARC36 DRIP36 VDRIP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	MED4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.