



KIF17 Antibody, FITC conjugated

Product Code	CSB-PA865208LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2E2
Immunogen	Recombinant Human Kinesin-like protein KIF17 protein (500-777AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF17 (KIF3-related motor protein), KIF17, KIAA1405 KIF3X
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF17
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.