



KIF17 Antibody, Biotin conjugated

Product Code	CSB-PA865208LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P2E2
Immunogen	Recombinant Human Kinesin-like protein KIF17 protein (500-777AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF17 (KIF3-related motor protein), KIF17, KIAA1405 KIF3X
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIF17