



POT1 Antibody, FITC conjugated

Product Code	CSB-PA868326LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUX5
Immunogen	Recombinant Human Protection of telomeres protein 1 protein (250-432AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protection of telomeres protein 1 (hPot1) (POT1-like telomere end-binding protein), POT1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	POT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.