



POLL Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA871513LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UGP5 |
| Immunogen | Recombinant Human DNA polymerase lambda protein (16-211AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | DNA polymerase lambda (Pol Lambda) (EC 2.7.7.7) (EC 4.2.99.-) (DNA polymerase beta-2) (Pol beta2) (DNA polymerase kappa), POLL |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | POLL |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |