




SLC25A23 Antibody

Product Code	CSB-PA871580LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV35
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-3 protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein SCaMC-3 (Mitochondrial ATP-Mg/Pi carrier protein 2) (Mitochondrial Ca(2+)-dependent solute carrier protein 2) (Small calcium-binding mitochondrial carrier protein 3) (Solute carrier family 25 member 23), SLC25A23, APC2 MCSC2 SCAMC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A23

Image



Western Blot

Positive WB detected in: Mouse brain tissue,
Mouse kidney tissue

All lanes: SLC25A23 antibody at 3µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 53, 50, 55, 49 kDa

Observed band size: 53 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.