



# SLC25A23 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA871580LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9BV35
<b>Immunogen</b>	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-3 protein (1-188AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Calcium-binding mitochondrial carrier protein SCaMC-3 (Mitochondrial ATP-Mg/Pi carrier protein 2) (Mitochondrial Ca(2+)-dependent solute carrier protein 2) (Small calcium-binding mitochondrial carrier protein 3) (Solute carrier family 25 member 23), SLC25A23, APC2 MCSC2 SCAMC3
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	SLC25A23
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.