



SLC25A23 Antibody, FITC conjugated

Product Code	CSB-PA871580LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV35
Immunogen	Recombinant Human Calcium-binding mitochondrial carrier protein SCaMC-3 protein (1-188AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium-binding mitochondrial carrier protein SCaMC-3 (Mitochondrial ATP-Mg/Pi carrier protein 2) (Mitochondrial Ca(2+)-dependent solute carrier protein 2) (Small calcium-binding mitochondrial carrier protein 3) (Solute carrier family 25 member 23), SLC25A23, APC2 MCSC2 SCAMC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC25A23