



PRDM16 Antibody, Biotin conjugated

Product Code	CSB-PA872534LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HAZ2
Immunogen	Recombinant Human PR domain zinc finger protein 16 protein (448-683AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PR domain zinc finger protein 16 (PR domain-containing protein 16) (Transcription factor MEL1) (MDS1/EVI1-like gene 1), PRDM16, KIAA1675 MEL1 PFM13
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PRDM16