



FOXQ1 Antibody, HRP conjugated

Product Code	CSB-PA874860LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C009
Immunogen	Recombinant Human Forkhead box protein Q1 protein (10-123AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein Q1 (HNF-3/forkhead-like protein 1) (HFH-1) (Hepatocyte nuclear factor 3 forkhead homolog 1), FOXQ1, HFH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXQ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.