



FOXQ1 Antibody, FITC conjugated

Product Code	CSB-PA874860LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C009
Immunogen	Recombinant Human Forkhead box protein Q1 protein (10-123AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Forkhead box protein Q1 (HNF-3/forkhead-like protein 1) (HFH-1) (Hepatocyte nuclear factor 3 forkhead homolog 1), FOXQ1, HFH1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FOXQ1