



GPAM Antibody, Biotin conjugated

Product Code	CSB-PA875724LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCL2
Immunogen	Recombinant Human Glycerol-3-phosphate acyltransferase 1, mitochondrial protein (281-459AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glycerol-3-phosphate acyltransferase 1, mitochondrial (GPAT-1) (EC 2.3.1.15), GPAM, GPAT1 KIAA1560
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPAM
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.